





# ORAL FLUID TEST


## RESULTS YOU TRUST



- **Compact.** Built-In Saturation Indicator Within Swab Collection Device
- **Secure.** Swab Locks In Place, Creating An Airtight Tamper Evident Seal
- **Configurable.** 16 Assay Strips Available For Configuration

RESULTS YOU CAN TRUST

POINT OF CARE MEDICAL  
DIAGNOSTIC TESTING PRODUCTS



# ORAL FLUID TEST

## RESULTS YOU TRUST

TEST	CALIBRATOR	CUT-OFF (NG/ML)
Amphetamine (AMP)	D-Amphetamine	50
Barbiturate (BAR)	Barbiturate	50/300
Benzodiazepine (BZO)	Oxazepam	10/50
Buprenorphine(BUP)	Buprenorphine	5/10
Cocaine (COC)	Cocaine	20
Cotinine (COT)	Cotinine	50
EDDP (EDDP)	2-Ethyliden-1,5-Dimethyl-3,3-Diphenylpyrrolidine	20
Methadone (MTD)	Methadone	30/75
Methamphetamine (MET)	D-Methamphetamine	50
Ecstasy (MDMA)	3,4-Methylenedioxymethamphetamine	50
Heroin 6-MAM	6-Monoacetylmorphine	25
Opiates (OPI)	Morphine	40
Oxycodone(OXY)	Oxycodone	40/20
Phencyclidine (PCP)	Phencyclidine	10
Propoxyphene(PPX)	Propoxyphene	50
Marijuana (THC)	11-nor-9-THC-9 COOH	12
Marijuana (THC)	9-THC	50
Alcohol (ALC)	Alcohol	>0.02"% B.A.C.

### ORAL FLUID TEST STOCKED CONFIGURATIONS

6 Panel  
THC/COC/OPI/mAMP/BZO/OXY

6 Panel  
THC/COC/OPI/mAMP/AMP/PCP

10 Panel  
THC/COC/OPI/mAMP/AMP/  
BZO/PCP/OXY/BAR/MTD